

### Product description

Order number	41.1395.105
Product description	Micro tube 1.3 ml, with screw cap, K3 EDTA, violet US code, paper label with lot & shelf life

### Product characteristics

Type of preparation	K3 EDTA
Preparation concentration	1.6 mg/ml blood
Colour code	US
Label/ Print	With paper label
Colour of print/label	Violet
Graduation	Yes
Type of cap	Screw cap
Base shape	Conical base with skirted base
Recommended centrifugation	2000 x g
Centrifugation time	10 min

### Size

Sample volume	1,3 ml
Diameter	10,8 mm
Length of product	47 mm
Length excluding cap	45 mm

### Material & colours

Product material	Polypropylene (PP)
Colour of product	Transparent
Cap material	Polypropylene (PP)
Colour of cap	Violet
Tube material	Polypropylene (PP)

### Purity & certification

Satisfies the requirement	IATA, ADR
Product in accordance with MDD or IVDD	IVD

CE certified	CE
Lot or serial number	batched

### **Packaging**

Minimum order qty.	1000
Pcs. / inner box	100
Pcs. / box	1000
Pcs. / pallet	96000
Depth of case	348 mm
Width of case	278 mm
Height of case	140 mm
Case volume	0,0135 cbm
Weight of case	1,89 KG
EAN of inner box	4038917128407
EAN of outer box	4038917021708